Flexi-Seal™ PROTECT PLUS

Luer Irrigation / Rx port

Fecal Management System / Stuhlmanagementsystem / Sistema di gestione dell'incontinenza fecale / Système de prise en charge de matières fécales / Sistema de control fecal / Sistema de controlo fecal / Fækalkatetersystem / Ulosteehallintajärjestelmä / Fekalkatetersystem / Fecaal managementsysteem / Система за управление на фекалиите / Sistem za upravljanje z blatom / Systém na odvádzanie stolice / Székletürítést szabályozó rendszer / Systém na odvod tekuté a polotekuté stolice / Väljjaheite kogumiskomplekt / Sustav za zbrinjavanje fecesa / System do kontrolowanej zbiórki stolca / Fëçu savākšanas sistēma / Fekalijū tvarkymo sistema / Sistemul de Management al fecalelor / Fekal Yönetim Sistemi / Σύστημα διαχείρισης κοπράνων / جهاز التحكم في البراز







Figure 1


Directions for Use / Gebrauchsanweisung / Istruzioni per l'uso / Mode d'emploi / Instrucciones de uso / Instruções de Utilização / Brugsanvisning / Käyttöohjeet / Bruksanvisning / Gebruiksaanwijzing / Uказания за употреба / Navodila za uporabo / Návod na použitie / Használati utasítás / Pokyny k použití / Kasutusjuhend / Upute za uporabu / Wskazówki dotyczące stosowania / Lietošanas norādījumi / Naudojimo instrukcija / Instrucțiuni de utilizare / Kullanım talimatları / Οδηγίες χρήσης / تعليمات الاستعمال

Fecal Management System Flexi-Seal™ PROTECT PLUS

Contents

(en) English	3	(sk) Slovenský	69
(de) Deutsch	8	(hu) Magyar	74
(it) Italiano	13	(cs) Český	79
(fr) Français	18	(et) Eesti	84
(es) Español	23	(hr) Hrvatski	89
(pt) Portuguesa	28	(pl) Polski	94
(da) Dansk	33	(lv) Latviešu	99
(fi) Suomi	38	(lt) Lietuvių	104
(no) Norsk	43	(ro) Română	109
(sv) Svenska	48	(tr) Türkçe	114
(nl) Nederlands	53	(el) Ελληνικός	119
(bg) Български	58	(ar) عربي	124
(sl) Slovenski	64		

- 1 Catheter Irrigation and Medication Port
- 2 Balloon Inflation Port with 2-Dome Fill Indicator and Cap
- 3 Sampling Port
- 4 Catheter
- 5 Position Indicator Line
- 6 Low-Pressure Retention Balloon with Finger Pocket
- 7 Self-Closing Connector
- 8 Hanging Strap
- 9 Luer-Lock Syringe
- 10 Collection Bag
- 11 Cinch Clamp
- 12 4 ConvaTec Diamonds™ Gelling and Odor Control Sachets

The Flexi-Seal™ PROTECT PLUS Fecal Management System contains (Figure 1):

- 1** privacy bag with filter
- 1** self-closing soft catheter tube assembly
- 1** Luer-Lock syringe
- 1** cinch clamp
- 4** ConvaTec Diamonds™ gelling and odor control sachets

The soft catheter is inserted into the rectum for fecal management to contain and divert fecal waste to protect the patient's skin and keep the bedding clean. There is a low-pressure retention balloon at the distal end and a connector for attaching the collection bag at the other end. There is a recess under the balloon for the clinician's finger allowing the device to be positioned digitally. A blue and a white port are attached to the side of the catheter. The white port, marked "≤45ml" is used to inflate the retention balloon after the device has been inserted into the patient's rectum. (Figure 2-1b). This white inflation port is equipped with two fill indicating domes (PROTECT indicator domes), green color (i.e., the dome closest to the catheter tubing) and red color (i.e., the dome furthest from the catheter tubing). The green fill indication dome provides a visual and tactile indication as to when the low-pressure retention balloon is filled to its optimal volume. A red-dome fill indicator provides a visual and tactile indication when the balloon over-inflation occurs. A white cap is provided to close off the white inflation port after the balloon inflation. The Luer-Lock connector in the blue housing, marked "IRRIG. /Rx" is used to flush the device if needed and administer medication, if prescribed. (Figure 2-1a). A dark blue sampling port is also present if stool samples are required to be taken by the clinician.

Intended Clinical Benefit

As determined by the available data, the clinical benefits of the Flexi-Seal™ FMS family of products include:

- Containment and diversion of faecal matter
- The enablement of being able to administer medication

Claims

In addition to the safety and performance claims relating to the safe use of the device, performance claims made for the device include:

- A product shelf life of 5 years
- The collection bag shall have 1 litre volume capacity.
- The Collection Bag shall have volume graduation marking ranging from 100 ml to 1,000 ml with an accuracy of ± 15%.

Indications for use

The Flexi-Seal™ PROTECT PLUS Fecal Management System is an indwelling fecal management catheter intended for use to manage fecal incontinence through the collection of liquid to semi-liquid stool and to provide access to administer medications. The device is intended for use in adult patients.

Contraindications

1. This product is not intended for use
 - for more than 29 consecutive days
 - for pediatric patients as its use has not been tested in this population
2. The Flexi-Seal™ PROTECT PLUS Fecal Management System should not be used on individuals who
 - have suspected or confirmed rectal mucosal impairment, i.e. severe proctitis, ischemic proctitis, mucosal ulcerations
 - have had rectal surgery within the last year
 - have any rectal or anal injury
 - have hemorrhoids of significant size and/or symptoms
 - have a rectal or anal stricture or stenosis
 - have a suspected or confirmed rectal/anal tumor
 - have any in-dwelling rectal or anal device (e.g. thermometer) or delivery mechanism (e.g. suppositories or enemas) in place
 - are sensitive to or who have had an allergic reaction to any component within the system

Warnings

- **Warning:** Clinicians should be aware that there are very limited clinical data on the use of in-dwelling fecal management systems after 14 days continued use.
- **Warning:** There is a potential risk of misconnections with connectors from other healthcare applications, such as intravenous equipment, breathing and driving gas systems, urethral/urinary, limb cuff inflation neuraxial devices and other enteral and gastric applications.
- **Warning:** Not following these instructions for use may increase the likelihood of an adverse event.
- **Warning:** Patients should be monitored daily and a physician notified immediately if any of the following occur
 - Rectal pain
 - Rectal bleeding
 - Abdominal symptoms such as distension/pain
- **Warning:** Over inflation of the retention balloon has the potential to increase the risk of adverse events including rectal pain, bleeding, ulcerations, and possible perforations.
- **Warning:** There is a danger of fecal obstructions with this product.
- **Warning:** The output may appear darker than usual and/or may contain black flecks. This is a visible indication of the ConvaTec Diamonds™. If monitoring output color, please use the sampling port or catheter. In case of contact with eyes, rinse immediately with clean water and seek medical advice. Store the ConvaTec Diamonds™ in a cool dry place. Do not open sachet. Keep out of the reach of children.

Precautions and observations

1. Close attention should be exercised with the use of the device in patients who have inflammatory bowel conditions or who have had rectal surgery. The physician should determine the degree and location of inflammation or extent of surgery (e.g. location of anastomosis) within the colon/rectum prior to considering use of this device in patients with such conditions.
2. Care should be exercised in using this device in patients who tend to bleed from either anti-coagulant / antiplatelet therapy or underlying disease. If signs of rectal bleeding occur, remove the device immediately and notify a physician.
3. The device should be used with caution in patients with spinal cord injury because of the possibility of the development of autonomic dysreflexia.
4. Remove any indwelling or anal device prior to insertion of the Flexi-Seal™ PROTECT PLUS and do not insert any other devices into the rectum while the Flexi-Seal™ PROTECT PLUS is in place.
5. Ensure that the patient does not lie or sit on the catheter as this could lead to localized pressure damage and contribute to the development of anal skin breakdown and/or restrict fecal flow.
6. Solid or soft-formed stool cannot pass through the catheter and will obstruct the opening. The use of the device is not indicated for solid or soft-formed stool.
7. Small amounts of moisture or seepage around the catheter is anticipated. To avoid skin irritation, initiate an appropriate institutional skin care protocol. At a minimum, the skin should be kept clean, dry and protected with a moisture barrier product.

8. If the catheter becomes blocked with feces, it can be rinsed with water using the irrigation port only (see Direction for Use "Irrigation of the Device"). DO NOT use the white inflation port (marked "≤45ml") to irrigate. If obstruction of the catheter is due to solid stool, use of the device should be discontinued.
9. Clinicians should take extra care to use the blue irrigation/medication housing (marked "IRRIG./Rx" Figure 7-1a) only when irrigating and delivering medication. DO NOT irrigate or administer medication through the white inflation port (marked "≤45ml" Figure 7-1b) which is closed off with a white cap.
10. Discontinue the use of the device if the patient's bowel control, consistency and frequency of stool begin to return to normal.
11. As long as the patient is regularly and closely monitored at all times, patients may be seated for short periods i.e. for up to 2 hours, as part of daily nursing care. During this period of seating, regular monitoring should be made to ensure the tubing is never blocked or kinked and to check for and avoid pressure damage to the anal/peri-anal region. For some patients, the length of the period of seating to avoid pressure damage to the anal/peri-anal region could be much shorter and the clinician should be alert to this possibility.
12. As with the use of any rectal device, the following adverse events could occur:
 - Leakage of stool around the device
 - Rectal/anal bleeding due to pressure necrosis or ulceration of rectal or anal mucosa
 - Peri-anal skin breakdown
 - Temporary loss of anal sphincter muscle tone
 - Infection
 - Bowel obstruction
 - Perforation of the bowel
13. This device is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.
14. If there is no flow of stools in a 24hr period then the following actions should occur, irrigation (see section D. irrigation of the device, in Directions for use) or removal of the device.
15. The collection bag should be coupled to the catheter in the correct orientation as shown in Fig. 3. Reading measurements of the collection bag is approximate only. The collection bag measuring scale has an accuracy of +/-15% and is for indication purposes only.
16. Do not use if package is damaged. Do not use Diamonds™ sachets if sachets are significantly broken.
17. If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority. "Circumstances in which the user should consult a health care professional"
18. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

Directions for use

A. Preparation of the device

1. In addition to the device system, gloves and lubricant will be required.
2. Unfold the length of the catheter to lay it flat on the bed, extending the collection bag toward the foot of the bed.
3. Insert 3 or 4 ConvaTec Diamonds™ sachets, one at a time, into the bag opening. Do not force the sachets, if resistance is met then gently move the sachet sideways ensuring they are placed at the bottom of the bag. (WARNING: DO NOT use the bag content for source of clinical information on stool color or consistency as it is modified by the gelling agent. Do not open sachet.)
4. Securely attach the collection bag to the connector at the end of the catheter with the correct orientation without trapping the bag against the bag connector (Figure 3). WARNING: trapping the bag against the bag connector would prevent coupling to the catheter.
 - 4.1 Follow the label shown on the catheter connector to couple the collection bag and catheter correctly.
 - 4.2 Position the catheter connector at a 90 degree angle to the bag connector opening and gently insert the catheter connector into the bag connector.
 - 4.3 Locate the two pins on the collection bag connector and align them to the two corresponding slots on the catheter connector.
 - 4.4 Gently push the catheter connector into the bag connector and then twist clockwise to securely attach the collection bag to the connector at the end of the catheter.

B. Preparation of the patient

1. Position the patient in left side lying position; if unable to tolerate, position the patient so access to the rectum

is possible.

2. Remove any in-dwelling or anal device prior to insertion of the Flexi-Seal™ PROTECT PLUS device.
3. Perform a digital rectal exam to evaluate suitability for insertion of device.

C. Insertion of the device and inflation of the balloon using a Luer syringe

1. Remove the white cap from the inflation port. Using the Luer syringe provided, remove the air that is in the balloon by attaching the Luer syringe to the white inflation port (marked “≤45ml”) and withdrawing the plunger. Remove the supplied Luer syringe and fill it with 45ml of water or saline and connect the Luer syringe to the white inflation port of the catheter. Insert a lubricated gloved finger into the blue finger pocket for digital guidance during device insertion (the finger pocket is located above the position indicator line). (Figure 4) Coat the balloon end of the catheter with lubricant. Grasp the catheter and gently insert the balloon end through the anal sphincter until the balloon is beyond the external orifice and well inside the rectal vault. The finger may be removed or remain in place in the rectum during initial balloon inflation.
2. Inflate the balloon with up to 45ml of fluid by slowly depressing the Luer syringe plunger. With the insertion finger removed, the green indication dome will indicate once the balloon has reached the optimal fill level for the anatomy. (Figure 5). Stop inflation once the green dome has signaled optimal fill. Under no circumstances should the balloon be inflated with more than 45ml of fluid. If the green indication dome indicates at less than 30ml of fluid, withdraw the fluid and reposition the balloon in the rectal vault. After repositioning, fill the balloon as described above. Do not fill with more than 45ml of fluid. If the red indication dome starts to inflate, assess the patient’s position, fully deflate the balloon and repeat the balloon inflation process. Stop inflation once green dome has signaled optimal fill.
3. Remove the Luer syringe from the inflation port, and gently pull on the soft catheter to check that the balloon is securely in the rectum and that it is positioned against the rectal floor. (Figure 6). Close the cap on the white inflation port to avoid misconnection issues.
4. Position the length of the flexible catheter along patient’s leg avoiding kinks and obstruction. Take note of the position indicator line relative to the patient’s anus. Regularly observe changes in the location of the position indicator line to determine movement of the retention balloon in the patient’s rectum. This may indicate the need for the balloon or device to be re-positioned. In the event of expulsion of the device, deflate the balloon fully; rinse the balloon end of the catheter and reinsert following the instructions for “Insertion of Device”. A rectal exam should be conducted prior to re-insertion to verify that no stool is present. If expulsion continues for more than three episodes, discontinuation of the device should be considered.
5. Hang the bag by the bead strap on the bedside at a position lower than that of the patient. Add insertion date of the fecal management system onto the supplied date formatted adhesive label and affix label onto flat region of the hanging strap.

D. Irrigation of the device

To irrigate the device, fill a Luer syringe with water at room temperature, attach the syringe to the connector in the blue irrigation/medication housing (marked “IRRIG./Rx” Figure 7-1a) and slowly depress the plunger. Do not irrigate through the white inflation port (marked “≤45ml” Figure 7-1b) as this would lead to over inflation of the retention balloon and the device would not be irrigated as intended. Repeat the irrigation procedure as often as necessary to maintain proper functioning of the device. Flushing the device as described above is an optional procedure for use only when needed to maintain the unobstructed flow of stool into the collection bag. If repeated flushing with water does not return the flow of stool through the catheter, the device should be inspected to ascertain that there is no external obstruction (i.e. pressure from a body part, piece of equipment, or resolution of diarrhea). If no source of obstruction of the device is detected, use of the device should be discontinued.

E. Maintenance of the device

Change the collection bag as needed per Step A.4. After removal of the bag from the catheter, close the bag with the provided cap. To remove the collection bag, push the catheter connector into the bag connector and then twist anti-clockwise to disengage. See correct method of holding the collection bag without trapping the bag against the bag connector in figure 3. Gently pull the catheter connector from the collection bag. Use the thumb to push the back of the cap into the bag connector which is to be held in place from the rear of the collection bag using the middle and index fingers. Use thumb to press around the cap to ensure full bag closure. Discard used bags per institutional protocol for disposal of medical waste. Observe the device frequently for obstructions from kinks, solid fecal particles or external pressure.

F. Administration of medication (dosage needs to be controlled by physician) using 2 syringes (1 for flushing with water, 1 for administering medication)

1. Fill a syringe with 10ml of water, attach the syringe to the blue connector and flush the irrigation line with 10ml of water.
2. Prepare a new syringe with prescribed medication. Position the cinch clamp loosely on the catheter at the black indicator line. Connect the syringe to the connector in the blue irrigation/medication housing ("IRRIG./R" Figure 7-1a) and administer medication. Clinicians should take extra care to use the connector in the blue irrigation/medication housing only when delivering medication. Do not administer medication through the white inflation port (marked " $\leq 45\text{ml}$ ") Figure 7-1b) as this would lead to over inflation of the retention balloon and the patient would not receive medication as intended.
3. To ensure delivery of medication into the rectum fill the syringe used for water flushing with 50ml of water and immediately flush the irrigation line. Tighten the cinch clamp on the catheter to ensure no flow through the catheter (ensure the second notch is engaged; squeeze tightly using forefinger and thumb of both hands to ensure a good seal). Allow the medication to dwell in the rectum for the desired amount of time as dictated by the prescribing physician.
4. Remove the cinch clamp. Fill the syringe used for water flushing with 10ml of water and flush the irrigation line. Dispose both syringes per institutional policy.

G. Sampling

To collect a sample from the catheter, open the dark blue sample port cap. Press the tip of a Luer-slip syringe or catheter syringe through the slit inside of the sampling port to access the interior of the catheter. Withdraw the syringe plunger to collect the sample. Withdraw the syringe and close the dark blue sampling port cap.

H. Removal of the device and deflation of the balloon

To remove the catheter from the rectum, the retention balloon must first be deflated. Remove the white cap from the inflation port. Attach a Luer syringe to the white inflation port (marked " $\leq 45\text{ml}$ ") and slowly withdraw all fluid from the retention balloon. (Figure 2-1b).

Disconnect the Luer syringe and discard. Grasp the catheter as close to the patient as possible and slowly remove from the anus. Dispose of the device in accordance with institutional protocol for disposal of medical waste. If the balloon is difficult or impossible to deflate, cut the inflation lumen and drain out the water in the balloon. On no account should the device be removed from a patient with the balloon still inflated.

General guidelines

The device may be changed as needed to perform normal patient assessment.

The device is not intended for use for more than 29 consecutive days.


MRI Safety Information

 Non-clinical testing has demonstrated that the Flexi-Seal™ PROTECT PLUS is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T or 3.0 T
 - Maximum spatial field gradient of 2,000 gauss/cm (20 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- The presence of this device may produce an image artifact.

©2023 ConvaTec Inc.

™ indicates trademarks of ConvaTec Inc. ConvaTec, the ConvaTec Logo, Flexi-Seal, the Flexi-Seal Logo are registered trademarks of ConvaTec Inc. in the U.S.

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CH REP

ConvaTec International Services GmbH
Mühentalstrasse 36/38
8200 Schaffhausen, Switzerland

ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

ConvaTec Turkey
Ayazağa Mah. Mimar Sinan Sk.
A Blok, No:21A, İç Kapı No: 9,
34396 Sarıyer
İstanbul TURKEY

Sponsored in Australia by:
ConvaTec (Australia) Pty Ltd.
PO Box 63, Mulgrave VIC 3170
Australia
1800 339 412 Australia
0800 441 763 New Zealand

© 2023 ConvaTec Inc.
® / ™ indicates trademarks of ConvaTec Inc.
ConvaTec, the ConvaTec Logo, Flexi-Seal and
the Flexi-Seal Logo are registered trademarks of
ConvaTec Inc. in the U.S.

Revised 2023-09
SI2352VA

Quick Reference for Device Insertion

ATTENTION: Please refer to the instructions for use for full information on device insertion.

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

- 1** Attach collection bag to catheter connector with arrows pointing upwards.


- 2** Connect Luer syringe to white inflation port and remove all air.


- 3** Lubricate finger and insert finger into blue finger pocket on balloon and lubricate balloon.

- 4** Gently insert the balloon into the rectal vault.

- 5** Connect the syringe to the white inflation port and inflate the balloon with water or saline until green indicator pops.

- 6** Gently pull on soft catheter to check secure placement against rectal floor.